

Artik Handheld Phototherapy

Prescription Models

OVERVIEW

The new handheld ultraviolet phototherapy devices are for prescription treatment of many photo-responsive skin disorders such as psoriasis and vitiligo. The devices are available with one or two Narrowband UVB lamps and includes an attachment that ensures a consistent distance to the skin.

KEY FEATURES

- Bright OLED controller display
- Automatic lamp shutoff at end of treatment
- FDA cleared
- Compatible with prescription protocols
- Covered by many insurance plans

LAMP(S)

- One or two Philips brand 9W 31 Inm dual-tube Narrowband UVB lamps (9 or 18 watts total UV energy)

INSURANCE

- May be eligible for insurance coverage subject to diagnosis and other factors
- HCPCS reimbursement code E0691
- Two lamp unit is an insurance upgrade



Includes comb/spacer, power cord and protective eye goggles



Rx only

INTERNATIONAL PRODUCT SPECS

- For countries with electrical current other than U.S. standard, a step-down transformer is required.

ELECTRICAL RATING

- 120 VAC, 60 Hz, 0.5A
- NEMA 1-15P polarized power cord

WARRANTY


- All parts and labor (excluding lamps) carry a 1-year warranty
- Lamps carry a 90-day warranty



Authorized Distributor
www.artikmed.com

Artik Medical Supply
888-433-9796
contact@artikmed.com

Prescription Handheld Models Comparison



**Single-Lamp Handheld
KN-4003BL2S**

DIMENSIONS

Output window: 1.6" x 5.1" = 8.2 in²
(4.1cm x 13.0cm = 53.3cm²)

Wand (no comb): 13" x 2.3" x 1.8"
(33.0cm x 5.7cm x 4.4cm)

Wand with Comb: 13" x 2.3" x 2.4"
(33.0cm x 5.7cm x 6.0cm)

AC adapter: 4.5" x 2.3" x 1.5"
(11.5cm x 5.9cm x 3.8cm)

WEIGHT

Wand only: 0.88 lbs. (0.4 kg)

Wand with power adapter and cord:
1.89 lbs. (0.86 kg)

APPLICATION

Good for a few small spots or distributed small spots. In some cases can be used on the scalp. Good for very curved skin surfaces such as arms and legs.



**Two-Lamp Handheld
KN-4006BL1S**

DIMENSIONS

Output window: 2.9" x 5.0" = 14.5 in²
(7.4cm x 12.7cm = 94cm²)

Wand (no spacer): 14.0" x 4.0" x 1.9"
(35.6cm x 10.2cm x 4.8cm)

Wand with spacer: 14.0" x 4.2" x 2.5"
(35.6cm x 10.7cm x 6.4cm)

AC adapter: 4.8" x 2.5" x 2.8"
(12.2cm x 6.4cm x 7.1cm)

WEIGHT

Wand only: 1.19 lbs. (0.54 kg)

Wand with power adapter and cord:
2.93 lbs. (1.3 kg)

APPLICATION

Can be used for small or distributed spots on curved skin, but also good for larger spots on flatter skin surfaces such as the back, chest, thigh, stomach, etc.